

## Experiment No.1

### Determination the ratio of water of hydration in crystalied $\text{BaCl}_2 \cdot \text{XH}_2\text{O}$

#### Introduction

The water crystallization can be determined using sample weight change as a result of its water content ejection by heating it.

This method can be only used in the following cases:

- 1- When water is the primary part to be volatilized only.
- 2- When the composition of the precipitate does not change as a result of Oxidation or hydrolysis.

The Barium chloride salt loses all crystallized water at a temperature range from (100-120) centigrade, high temperatures can be used for this purpose considering that the salt does not break at high temperatures up to (800-900) centigrade.



Barium dichloride crystals shows a little tendency to absorb water in the spaces of the molecules of the substance, which means that any change in weight represents a great precision water of crystallization.

Crystallization water can be determinte for each of the following salts:



## Materials and tools

- 1- Crucible.
- 2- Desiccator.
- 3- Oven & Furnace.
- 4- Aqueous Barium Chloride crystals  $\text{BaCl}_2 \cdot X\text{H}_2\text{O}$  .

## Procedure

- 1- Heat the crucible inside the oven and then leave for a (30) minutes , cool in a desiccator for (15) minutes and then weigh after cooling.
- 2- Put in the crucible (0.5) gm of the sample and weigh it again ( aqueous Barium Chloride  $\text{BaCl}_2 \cdot X\text{H}_2\text{O}$  ).
- 3- Put the crucible inside the burning furnace for half an hour.
- 4- Let the crucible and its contents to be cooled in a desiccator and then weight it.
- 5- From weight loss after burning process, calculate the prcentage of water in aqueous Barium Chloride and Calculate the water molecule number in aqueous Barium Chloride ( $\text{BaCl}_2 \cdot X\text{H}_2\text{O}$ ) .

## Calculation

Weight of Water salt (model weight) = 0.5 g

Weight of anhydrous Salt = Weight of crucible with anhydrous Salt

(after heating) - Weight of an empty Crucible

The weight of water = the weight of aqueous salt - the weight of anhydrous salt

$$100 \times \frac{\text{wt of H}_2\text{O}}{\text{wt of sample}} = \text{H}_2\text{O}\%$$

Calculate the number of water molecules crystallization (x)

It is extracted from the following equation:

$$\frac{\text{M. wtBaCl}_2 \cdot \text{XH}_2\text{O}}{\text{M. wtBaCl}_2} = \frac{\text{wt BaCl}_2 \cdot \text{XH}_2\text{O}}{\text{wt BaCl}_2}$$

## Questions

- 1- What are the conditions required in a materials that water crystallization can be determined? Give an example.
- 2- Why empty Crucibles need to be heated?
- 3- Define crystallized materials? Can some of them be used in chemical analyzes?
- 4- What is the purpose of the crucible cooling after heating in a desiccator?