

## **Experiment (7)**

### **Application of (Lambert – Beer) law**

#### **Introduction:**

Spectrophotometric methods include measuring the intensity of incident light and transmitted light at a specific wavelength. There are two types of device are the spectral:  
(1) Single beam (2) Double beam.

#### **Law (Lambert – Beer):**

Absorbs equal parts of the light beams by equal changes in the concentration of the absorbed material when the optical path length in the absorbing material constant.

$$A = \epsilon C l$$

$$\%T = (I / I_o) \times 100$$

$$A = \log I_o / I$$

$$= \log 1 / T$$

$A$  = absorbance

$\epsilon$  = molar absorbance coefficient ( $\text{mol}^{-1} \cdot \text{L} \cdot \text{cm}^{-1}$ )

$l$  = long the light path inside the solution (1cm)

$T$  = transmittance

$C$  = concentration ( $\text{mol} \cdot \text{L}^{-1}$ )

$I$  = the intensity of transmitted light

$I_o$  = intensity of incident light

**Wavelength ( $\lambda$ ):** The distance between adjacent peaks of a wave packet, that has the standard unit:  $\text{nm} = 10^{-9}\text{m}$ ,  $\text{\AA} = 10^{-10}\text{m}$ ,  $\mu\text{m} = 10^{-6}\text{m}$

**The greatest wavelength ( $\lambda_{\text{max}}$ ):** The wavelength, which has the highest absorption of the substance.

**Light intensity (I):** The number of photons absorbed per second.

**Photon:** Is units of energy.

#### **Beer-Lambert Law Applications**

This law finds applications in various fields such as:

- **Analytical chemistry:** This analysis mainly concentrates on the separation, quantification and identification of matter by spectrophotometry. There is no involvement of extensive pre-processing of the sample to get the results. For example, bilirubin count in a blood sample can be determined by using a spectrophotometer.

- **Atomic Absorption:** The application of atomic absorption spectrometry (AAS) for the determination of metal concentrations.
- **In atmosphere:** Solar or stellar radiation in the atmosphere can be described using this law. The law in atmospheric applications has a modified equation.


### **Beer-Lambert Law Limitations**

Using this law it becomes easy to study the absorptivity coefficient of the sample when the concentration is low  $<10\text{mM}$  but as the concentration becomes high  $>10\text{mM}$  there is a deviation as the electrostatic interactions become more. The value of absorbance should be between (0-1).


### ***Measuring the absorbance of potassium permanganate ( $\text{KMnO}_4$ )***

#### **a) Finding ( $\lambda_{\text{max}}$ ):**

- 1) Prepare dilute solutions of  $\text{KMnO}_4$  (0.1M):  $(1, 3, 5, 7) \times 10^{-4}\text{M}$  in volumetric flask (50ml) according to law ( $M_1 V_1 = M_2 V_2$ ).
- 2) Measure absorbance of the lower concentration ( $1 \times 10^{-4}\text{M}$ ) versus water (blank) in the wavelength range (400 – 600 )nm and that each (10nm) to note the highest absorption value read by the device and install it to be ( $\lambda_{\text{max}}$ ).




Wave length(nm)	Absorbance
400	
410	
-----	
-----	
600	



**b) Finding ( $\epsilon$ ) & Concentration of unknown:**

- 1- When you install ( $\lambda_{max}$ ) from step (a) was appointed absorbance of each solution.
- 2- Plot a relationship between absorbance (A) and concentration then determine the molar absorption coefficient ( $\epsilon$ ) from the slope and the concentration of unknown.

Conc.(mol/L)	Absorbance
$1 \times 10^{-4}$	
$3 \times 10^{-4}$	
$5 \times 10^{-4}$	
$7 \times 10^{-4}$	



**Discussion:**

- 1- Is it possible to measure the potassium permanganate in the ultraviolet region near the visible region?
- 2- What is the difference between ultraviolet and visible spectrum?