Al-Karkh University of Science College of science First year level General chemistry Labs

Supervisor. Assis. Prof. Dr. Mohammed A. B. Abdul Jabar

Lab -9-

Separation of I Group cations (Ag⁺, Pb²⁺, Hg₂²⁺)

1. Introduction:

Group I cations (Ag⁺, Hg₂²⁺, Pb²⁺) form insoluble chlorides. Upon the addition of hydrochloric acid Ag⁺, Pb²⁺, Hg₂²⁺ ions will precipitate as AgCl, PbCl₂ and Hg₂Cl₂

AgNO ₃ + HCl → AgCl _↓ + HNO ₃	(reaction 1)	
$Pb(NO_3)_2 + 2 HCl \rightarrow PbCl_2\downarrow + 2 HNO_3$	(reaction 2)	
$Hg_2(NO_3)_2 + 2 HCI \rightarrow Hg_2Cl_2\downarrow + 2 HNO_3$	(reaction 3)	

The solubility of PbCl₂ increases approximately threefold as the temperature of the solution increases from 20°C to 100°C. Thus, PbCl₂ will dissolve in hot water while AgCl and Hg_2Cl_2 remains insoluble. The presence of Pb²⁺ ions in obtained solution can be proved with Kl, K_2CrO_4 and H_2SO_4 solutions. Their addition yield a golden yellow precipitates of Pbl₂ and PbCrO₄ (reactions 4 and 5), and white precipitate of PbSO₄ (reaction 6).

$Pb^{2+} + 2 KI \rightarrow Pbl_2 \downarrow + 2 K^{+}$	(reaction 4)
$Pb^{2+} + K_2CrO_4 \rightarrow PbCrO_{4\downarrow} + 2K^+$	(reaction 5)
$Pb^{2+} + H_2SO_4 \rightarrow PbSO_4 + 2H^+$	(reaction 6)

The precipitate still may contain Hg₂Cl₂ and AgCl. Of those two compounds, only the silver chloride is soluble in aqueous ammonia due to the formation a colorless solution of Ag(NH₃)₂Cl (reaction 7), whereas mercury(I) chloride turns into Hg metal and HgNH₂Cl visible as grayish-black precipitate which is insoluble in ammonia solution (reaction 8)

The formation of white precipitate of AgCl in reaction of Ag(NH₃)₂Cl with diluted HNO₃ proves the presence of silver (reaction 9). The additional confirmation of silver ions is yellowish precipitate of silver iodide Agl obtained in the reaction with potassium iodide Kl (reaction 10).

$$Ag(NH_3)_2CI + 2 HNO_3 \rightarrow AgCI_{\downarrow} + 2 NH_4NO_3$$
 (reaction 9)
 $Ag(NH_3)_2CI + KI \rightarrow AgI_{\downarrow} + 2 NH_3 + KCI$ (reaction 10)

Scanned with CamScanner

2. Procedure.

Take 2 ml of unknown for identification and add drops of 6M HCl. You should obtain white precipitate (Precipitate A).

Transfer the precipitate to a small beaker, and boil with 1-5 ml water. Filter hot.

Wash the ppt. with hot wat Add NH ₄ OH to the precip		
Residue: black precipitate HgNH ₂ Cl + Hg	Precipitate A	Cool a portion of the solution: a white crystalline ppt. of PbCl ₂ is obtained if Pb is present in any quantity.
Hg ₂ ²⁺	Ag ⁺	Pb ²⁺
Filtrace 2 Pacts P	Divide the filtrate into two parts: a.) Acidify with dilute HNO3. White ppt. of AgCl. b.) Add a few drops of KI solution. Pale yellow ppt. of AgI. Ag ⁺ present.	Divide the hot filtrate into three parts: a.) Add K ₂ CrO ₄ solution. Yellow ppt. of PbCrO ₄ , insoluble in dilute acetic acid b.) Add KI solution. Yellow ppt. of PbI ₂ , soluble in hot water to a colourless solution which deposites brilliant yellow crystals upon cooling. c.) Add dilute H ₂ SO ₄ . White ppt. of PbSO ₄ , soluble in ammonium tartrate solution. Pb ²⁺ present.

3. Results:

4. Discuss:

- a- PbCl₂ soluble in hot water.
- b- Presence of Silver ion (Ag+) in the unknown solution.
- c- Divide the colorless filtrate solution of Ag(NH₃)₂Cl into two parts.

Take 2 ml of unknown for identification and add drops of 6M HCl. You should obtain white precipitate (Precipitate A).


Transfer the precipitate to a small beaker, and boil with 1-5 ml water. Filter hot.

Residue: AgCl, Hg ₂ Cl ₂ Wash the ppt. with hot wate Add NH ₄ OH to the precipi		Filtrate: PbCl ₂
Residue: black precipitate HgNH ₂ Cl + Hg	Solution: Ag ⁺ complex [Ag(NH ₃) ₂] ⁺	Cool a portion of the solution: a white crystalline ppt. of PbCl ₂ is obtained if Pb is present in any quantity.
Hg2 ²⁺	Ag ⁺	Pb ²⁺
2 1209	Divide the filtrate into two parts: a.) Acidify with dilute HNO3. White ppt. of AgCl. b.) Add a few drops of KI solution. Pale yellow ppt. of AgI. Ag+ present.	Divide the hot filtrate into three parts: a.) Add K ₂ CrO ₄ solution. Yellow ppt. of PbCrO ₄ , insoluble in dilute acetic acid b.) Add KI solution. Yellow ppt. of PbI ₂ , soluble in hot water to a colourless solution which deposites brilliant yellow crystals upon cooling. c.) Add dilute H ₂ SO ₄ . White ppt. of PbSO ₄ , soluble in ammonium tartrate solution. Pb ²⁺ present.

3. Results:

4. Discuss:

- a- PbCl2 soluble in hot water.
- b- Presence of Silver ion (Ag+) in the unknown solution.
- c- Divide the colorless filtrate solution of Ag(NH₃)₂Cl into two parts.



Scheme 1. Analysis of the Group I.